



BODILY FLUIDS TEST KITS LABORATORY + CRIME SCENE

SERAQUANT READER

WWW.SERATEC.COM

Store at 2-30°C

ABOUT US

Founded in 1987 as a research spin-off company in the old university town of Goettingen, Germany, SERATEC® GmbH has more than 25 years of experience in research and development, production and sales of immunochromatographic assays with particular focus on Forensic Biology.

Our test kits for the rapid detection of human bodily fluids help preselect circumstantial evidence for DNA analysis. Various products for both laboratory and crime scene testing are available to meet the sophisticated requirements of crime laboratories. Whatever we do, we make highest demands on ourselves. SERATEC® is ISO 9001 and ISO 13485 certified. Our employees face challenges of R&D, production and customer service with utmost professionalism.

Customer's opinion is the most valuable asset. We encourage forensic scientists from crime laboratories and educational institutions to be actively involved in product development and validation.

Consistently high quality, comprehensive knowledge and expertise, fair prices and great respect towards customers, partners and competitors. This is what SERATEC® stands for.





BODILY FLUIDS TEST KITS Today's forensic biologist needs a rapid, accurate, sensitive and inexpensive screening tool to analyze potential biological fluid stains and send them for DNA analysis. SERATEC® offers test kits for the detection of human seminal fluid, blood and saliva. SERATEC SERATEC

HUMAN SEMINAL FLUID



PSA Semiquant LABORATORY KIT

Cat.No. PSM400F/40

UNIT

Box, 40 cards, pouched with Pasteur pipette and desiccant in each pouch;
1 vial 50 mL of dilution buffer

- + Ready to use, robust kit; both cards and buffer can be stored at room temperature
- + Long shelf life of 24 months after manufacture Fast and accurate: reliable results after only 5-10 minutes
- + High sensitivity: detects 1 ng/mL PSA or a 10⁻⁶ dilution of semen
- + High specificity: cross reacts only with the seminal fluid of primates

PSA Semiquant CRIME SCENE KIT

Cat.No. PSM-CS/20

BOX 20 kits of 1 card (pouched, with
Pasteur pipette and desiccant)
1 vial with buffer solution
1 DNA-free swab for sample collection





HemDirect CRIME SCENE KIT

Cat.No. HBF-CS/20

BOX 20 kits of 1 card (pouched, with Pasteur pipette and desiccant)
1 vial with buffer solution

1 DNA-free swab for sample collection

+ uses 2 murine monoclonal anti-human hemoglobin antibodies as active components. When hemoglobin is present, a red line will appear in the test region for a positive result.

Internal control results is a red line in the

control region for a valid test.

- + Ready to use, robust kit; both cards and buffer can be stored at room temperature
- + Long shelf life of 24 months after manufacture
 Fast and accurate: reliable results after only
 5–10 minutes
- + Very high sensitivity: detects 20 ng/mL hemoglobin or a 10-7 dilution of blood
- High specificity: cross-reactivity only with upper primate blood and ferret blood (weak);
 no reaction with animal blood



AmylaseTest LABORATORY KIT

Cat.No. AMY/40

UNIT

Box, 40 cards, pouched with Pasteur pipette and desiccant in each pouch; 1 vial 50 mL of dilution buffer

- + Ready to use, robust kit; both cards and buffer can be stored at room temperature
- + Long shelf life of 24 months after manufacture
- + Fast and accurate: reliable results after only 5-10 minutes
- + Very Sensitive: 1:10,000 dilution of saliva shows positive result
- + Specific: no reaction with animal amylase except for primates; cross-reaction observed with pan creatic amylase

CRIME SCENE KIT

Cat.No. AMY-CS/20

BOX 20 kits of 1 card (pouched, with Pasteur pipette and desiccant) 1 vial with buffer solution

1 DNA-free swab for sample collection



THE SERAQUANT READER



The SeraQuant is an instrument for the computer-based, quantitative analysis of the SERATEC® test kits and the automatic, consistent result documentation.

The color intensity of the cassette result lines varies depending on concentration which enables a calculation of the concentration value (measurement) based on calibration curve.

With the SeraQuant, it is possible to create calibration curves for all three parameters used in Forensic Biology (PSA, Hemoglobin and alpha-Amylase) and therefore set user-specific standards of calibration and measurement quality.

The results of the cassette testing are automatically read by a camera under standard and consistent conditions and not by the naked eye. This procedure ensures that there is no room for subjective interpretation of results.

A unique sophisticated algorithm is able to accurately distinguish between very low positive and negative results.

Measurements are stored as secure files along with both image and descriptive data provided by user as well as calibration details for full traceability.

