SERATEC α -Amylase Test

Cat No: AMY
Rev. 08.2013

DE	Bedienungsanleitung
EN	Instructions for use
ES	Instrucciones de uso
FR	Mode d'emploi
IT	Istruzioni per l'uso
PL	Instrukcja użytkowania
RU	Инструкция по применению

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
Instructions for use: SERATEC α -Amylase Test Cat No: AMY

Intended Use

The SERATEC® Saliva test is a chromato-graphic immuno-assay for the rapid detection of human α -Amylase in forensic samples. It contains two monoclonal murine anti-human- α -Amylase antibodies as active compounds.

Introduction

Although found in many tissues, amylase is most prominent in saliva and pancreatic juice. There are at least four variants of human alpha-amylase, two of which are found in saliva, and the other two are secreted in the pancreas. Salivary α -Amylase breaks starch into maltose and dextrin. This form of amylase is also called ptyalin. It will break large, insoluble starch molecules into soluble starches (amylodextrin, erythrodextrin, and achrodextrin), producing successively smaller starches and ultimately maltose.

Presumptive test for saliva that detects the enzymatic activity of alpha-amylase (e.g. Phadebas), and not the alpha-amylase molecule itself, will yield a positive result if any alpha-amylase enzyme is present, regardless of the organism it came from.

The SERATEC® Saliva test eliminates these problems and shows a high specificity and sensitivity. Accordingly, the product is suitable for the identification of human saliva traces in forensic samples. It has the following advantages over conventional methods:

- Easy in use, directly at crime scenes or in the laboratory
- Quick and reliable result after only 10 minutes
- High sensitivity. Samples containing only 80 mIU/ml human α -Amylase react positive. Fresh human saliva diluted in the extraction buffer up to 1/2000 shows clear positive test results.
- The test is specific to human α -Amylase. Cross reactivity may occur with primate saliva.

α -Amylase may also be present in other body fluids (e.g. urine, stool, seminal fluid, vaginal fluid or blood), which may lead to positive results. However, the proper usage of SERATEC® Saliva test requires obligatory dilution of collected sample material. In this case, the test react positive with saliva samples only, whereas the α -Amylase concentration of other body fluids decreases below the detection limit.

A High dose Hook effect (prozone Effect) is not observed. Even neat saliva react positive with test, although it is not recommended to use neat saliva due to its high viscosity.

Materials provided

- 40 individually sealed α -Amylase cassettes, with one plastic pipette each
- 50 mL of standard buffer solution
- user instruction leaflet

Materials required but not provided: Timer

Storage and Stability

Both cassettes and buffer are stable up to the batch expiration dates stated on the sealed pouch and the buffer tube. Cassettes and buffer can be stored at room temperature or refrigerated (2-30°C resp. 38-86°F). The cassette must remain in the sealed pouch until use.

Qualitative Characteristics

Sensitivity: result line is visible at as low as 80 mIU/mL α -Amylase

The High dose Hook effect is not observed: human saliva is tested positive in the dilution range of 1:1 to 1:10³ using the provided buffer.

Reference Preparations

The qualitative characteristics of the test are confirmed in a final QC testing using the following standard: Sigma Aldrich α -Amylase from human saliva.

Specificity

The product shows no cross reactivity with other salivary proteins. No cross reactivity was observed with the saliva of the following animals: dog, cat, horse, cow, domestic pig, mouse, goat, hamster and rabbit.

It is verisimilar that saliva from upper primates reacts positive with the test.

Body fluids (blood, urine, seminal fluid) collected with a swab and solved in the provided buffer are reacting negative with the test.

Precautions

Forensic samples and all materials coming in contact with it should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire. The product itself and all materials coming in contact with forensic samples should be autoclaved before disposal.

- Do not use cassettes or buffer after expiration date.
- Do not use cassettes if the pouch has been damaged.
- The product consists of potentially infectious materials (e.g. antibodies). These materials do not cause any danger if the product is used according to the instructions.
- Do not open pouches until ready to perform the assay.
- Do not freeze any components of this product.

General Hints

- It is strongly recommended to use the provided buffer. It is adapted for obtaining optimal results. Other buffers may lead to varying line intensities.
- Do not use liquids with a pH-value below 3 or above 12 for testing. It can cause false or invalid results.
- Tissue particles do not interfere with the test result.
- Samples being subject to storage of one day or longer should be kept dry and cold (2-8°C). Liquid samples could be frozen.
- Viscous samples should be centrifuged for about one minute at 13,000g. The supernatant is used for testing. If the supernatant is still too viscous it should be further diluted.
- Cotton swabs and fabric cuttings should be extracted in an adequate amount of the provided dilution buffer. E.g., a cotton swab or a 1 cm² cutting could be extracted in 1mL of buffer.
- The addition of neat liquid samples from any source is not recommended. Liquid samples should be diluted at least 1:100 prior to testing.

Instructions for use: SERATEC α -Amylase Test Cat No: AMY


Start of the assay

- Allow all test components to warm up to room temperature before starting the test.
- Remove the cassette out of the foil pouch and tag the cassette with a marker, if necessary.
- Add three drops (about 120 μ l) in the sample well. Keep remaining sample if possible, in case it might be necessary to test additional dilutions.
- Read result after 10 minutes incubation time at room temperature. There should be no remaining fluid in the sample well at this time point.

Suggested reading/References

1. Evaluation of the SERATEC α -Amylase Saliva Test

Interpretation of results

Negative samples will show one line, the Control line (C), whereas α -Amylase positive samples will show two lines:

Test result line (T): reflects α -Amylase concentration of the sample, visible in α -Amylase-positive samples only

Control Line (C): control for possible procedural errors and for the integrity of test components

Negative result (no α -Amylase in the sample or α -Amylase concentration below detection limit)

Test result line (T) is not visible. Appearance of the control line (C) confirms validity of the test. In this case the sample most likely does not contain α -Amylase.

Note:

Make sure that the dilution of the probe leads to a α -Amylase concentration within the detection range. α -Amylase concentrations that are too low (e.g. due to insufficient extraction) may lead to false results.


Positive result (α -Amylase detectable)

Test result line (T) and control line (C) appear.

In this case it is very likely that the sample contains α -Amylase.

Invalid result

Control line (C) is not visible. The test is invalid and the assay should be repeated.

Symbols

	Expiration date
	Storage temperatur
LOT	Lot number

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