SERATEC HemDirect

Cat No: HBF07

Rev. 08.2013

DE	Bedienungsanleitung
EN	Instructions for use
ES	Instrucciones de uso
FR	Mode d'Emploi
IT	Istruzioni per l'uso
PL	Instrukcja użytkowania
RU	Инструкция по применению

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
Literatur

1. Amanda Misencik and Dale L. Laux: Validation Study of the Seratec HemDirect Hemoglobin Assay for the Forensic Identification of Human Blood. Midwestern Association of Forensic Science, Newsletter Spring 2007, p. 25 et seq.
2. Hochmeister et al. (1999) Validation Studies of an Immunochromatographic 1-Step Test for the Forensic Identification of Human Blood. J Forensic Sci Vol 44: 597-601.

Symbole

	Verfallsdatum
	Lagertemperatur
	Chargennummer

Rev: 08.2013

Intended use

The SERATEC® HemDirect Hemoglobin test is a chromatographic immunoassay for the rapid detection of human hemoglobin in forensic samples. It contains two monoclonal murine anti-human-hemoglobin antibodies as active compounds.

Introduction

Hemoglobin (Hb) is located in the erythrocytes and predominantly serves to carry oxygen and carbon dioxide within the body. Hemoglobin has a molecular weight of 64.5 kDa and is composed of 4 amino acid chains, two of which are identical. Each chain is associated with a hem group that contributes to around 4% of the total weight of the molecule. With concentrations between 120-160 mg/ml and 140-180 mg/ml for women and men, respectively, hemoglobin is one of the most abundant proteins in blood. Various confirmatory tests for the identification of blood are based on the detection of hemoglobin but do not allow differentiating between human and non-human blood. This question requires additional testing.

The SERATEC® HemDirect test eliminates these problems and combines a high specificity with an excellent sensitivity. Accordingly, the product is suitable for the identification of human blood traces in forensic samples. It has the following advantages over conventional methods:

- Ready-to-use product kits, no further equipment needed
- Easy in use, directly at crime scenes or in the laboratory
- Quick and reliable result after only 5-10 minutes
- Extremely high sensitivity. Samples containing only 40 ng/ml human hemoglobin react positive. Fresh human blood diluted in the extraction buffer up to 10^{-7} shows clear positive test results.
- The test is specific to human hemoglobin. Cross reactivity occurs with primate and ferret blood.
- Due to the long stability of hemoglobin, SERATEC® HemDirect is able to confirm presence of blood in very old stains. A study by Laux et al. demonstrates that the SERATEC® HemDirect is suitable for detecting hemoglobin in aged blood stains. Even a 31 year old stain reacted positive.¹

Dried blood stains and muscle tissue stored at different environmental conditions for one month reacted positive upon extraction. In contrast, additions of detergents like SDS, sarcosyl or bleach to the sample led to negative or weak positive results, probably caused by the denaturation of hemoglobin.

SERATEC® HemDirect results will be affected in their evidence by two facts. First, hemoglobin may also be present in other body fluids (e.g. urine, stool, seminal fluid, vaginal fluid or saliva), which may lead to positive results. However, the proper usage of SERATEC® HemDirect requires obligatory dilution of collected sample material. In this case, the test react positive with blood samples only, whereas the hemoglobin concentration of other body fluids decreases below the detection limit.

Second, very high concentration of human hemoglobin in the sample may lead to negative test results – a phenomenon known as high dose hook effect (prozone effect). This effect can be avoided by using the buffer and if necessary testing different dilutions of one sample. As the visible color caused by the hemoglobin vanishes between a 10^{-3} and a 10^{-4} dilution, the color of the sample might be a guide for the user. At these concentrations no high dose hook effect is observed.

Materials provided

- 30 individually sealed HemDirect cassettes, with one plastic pipette each
- 30 buffer tubes containing 1.5 mL of standard buffer solution
- User instruction leaflet

Materials required but not provided: Timer

Storage and Stability

Both cassettes and buffer are stable up to the batch expiration dates stated on the sealed pouch and the buffer tube. Cassettes and buffer can be stored at room temperature or refrigerated (+2-+30°C resp. 38-86°F). The cassette must remain in the sealed pouch until use.

Qualitative Characteristics

Sensitivity: result line is visible at as low as 40ng/mL Hem.

Fresh bloodstains should be diluted prior to testing to avoid false negative results caused by the high dose hook effect. Using the standard buffer, dilutions in the range from 1:50 to $1:10^7$ are tested positive. For this reason the recommended minimum dilution for fresh blood samples is 1:500.

Reference Material

The qualitative characteristics of the test are confirmed in a final QC testing using human hemoglobin purchased at SIGMA. The extraction buffer is used as solvent.

Specificity

The SERATEC® HemDirect Test shows no cross reactivity with bovine, dog, rabbit, cat, pig, wild boar, horse, chicken, sheep, mule, goat and red deer hemoglobin. Primate and ferret blood may react positive. Ref.1 shows the complete list.

Precautions

Forensic samples and all materials coming in contact with it should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire. The product itself and all materials coming in contact with forensic samples should be autoclaved before disposal.

- Do not use cassettes or buffer after expiration date.
- Do not use cassettes if the pouch has been damaged.
- The product consists of potentially infectious materials (e.g. antibodies). These materials do not cause any danger if the product is used according to the instructions.
- Do not open pouches until ready to perform the assay.
- Do not freeze the test cassettes.

Specimen collection and handling

- It is strongly recommended to use the provided buffer. It is adapted for obtaining optimal results. Other buffers may lead to varying line intensities and e.g. water to a decrease of sensitivity.
- Do not use liquids with a pH-value below 4 or above 12 for testing. It can cause false or invalid results. The application of the provided buffer will adjust the sample to an optimal pH-value.
- Tissue particles do not interfere with the test result.
- Samples being subject to storage of two days or longer should be kept dry and cold (2-8°C). Liquid samples could be frozen.
- Very old or small stains should be extracted longer than fresh and strong stains.¹
- Pure blood samples and samples of very high blood concentration should be diluted at least 1:500 prior to use.
- Viscous samples should be centrifuged for about one minute at 13,000g. The supernatant is used for testing. If the supernatant is still too viscous it should be further diluted.

Instructions for use: SERATEC HemDirect Cat No: HBF07

- Fabric cuttings from dried stains should be between 0.25-1 cm² for extraction with the buffer.
- The addition of neat liquid samples from any source is not recommended. Liquid samples should be diluted at least 1:500 prior to testing.
- No detergents like SDS or Sarcosyl should be added to the extraction buffer.

Test procedure

Specimen preparation

Hold collection tube upright and unscrew the purple cap. Collect the sample material (e.g. dried blood, blood stained thread of clothing) with the applicator stick attached at the cap or with other appropriate tools. Samples (e.g. cotton swabs, fabric cuttings) can be placed directly in the tube. Screw the cap tightly and extract the blood by shaking thoroughly. At this stage the specimen is stable at room temperature for at least 2 days.

The extraction of old bloodstains may be difficult.¹ If necessary the extraction can be done in the laboratory with the help of a shaker.


As solvent we recommend using the extraction buffer provided in the kit. However, other buffers at a neutral pH range (e.g. Phosphate, HEPES or TRIS buffered saline) might be used as well.

Extraction Buffer

The extraction buffer is composed of the following ingredients: 12.1 g TRIS; 8.8 g Na₃Citrate; 0.2 g NaN₃ → adjust pH auf 6.8 with HCl, 0.5 g Tween 20; 5.0 g BSA → fill with distilled water to 1000 ml

Assay Procedure

- Allow all test components to warm up to room temperature before starting the test. Only this will guarantee the stated sensitivity. A sample temperature of 8 °C leads to a decrease of the sensitivity to approximately 200 ng/ml.
- Remove the cassette out of the foil pouch and tag the cassette with a marker, if necessary.
- Add about 120 µl (three drops with the provided pipette) into the sample well of the cassette. Either directly out the buffer tube after breaking the seal of the tube by the use of a piece of paper (danger of splattering) or via a pipette. Keep remaining sample if possible, in case it might be necessary to test additional dilutions.
- Read result after 5 minutes incubation time at room temperature.
- Negative results should be confirmed after 10 minutes.

Interpretation of results

Negative result (no hemoglobin in the sample or hemoglobin below the detection limit)

Only one red colored line appears in the control region (C) after 10 minutes. The absence of a line in the test result region (T) indicates a negative test result.

Note:

Make sure that the dilution of the sample is within the detectable range of Hemoglobin concentrations!


Positive result (human hemoglobin present in the sample)

Two red colored lines appear in the result window, one at the control region (C) and one in the test result region (T). The color intensities of the two lines may vary. Even a weak test result line after 10 minutes indicates a positive test result.

Invalid result

The control line (C) does not appear.

In this case the test is invalid and should be repeated.

Literature

1. Amanda Misencik and Dale L. Laux: Validation Study of the Seratec HemDirect Hemoglobin Assay for the Forensic Identification of Human Blood. Midwestern Association of Forensic Science, Newsletter Spring 2007, p. 25 et seq.
2. Hochmeister et al. (1999) Validation Studies of an Immunochromatographic 1-Step Test for the Forensic Identification of Human Blood. J Forensic Sci Vol 44: 597-601.

Symbols

	Expiration date
	Storage temperatur
LOT	Lot number

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