# SERATEC PSA Semiquant

Cat No: PSM400F

Rev. 08.2013

<b>DE</b>	Bedienungsanleitung
<b>EN</b>	Instructions for use
<b>ES</b>	Instrucciones de uso
<b>FR</b>	Mode d'Emploi
<b>IT</b>	Istruzioni per l'uso
<b>PL</b>	Instrukcja użytkowania
<b>RU</b>	Инструкция по применению

## Anwendung

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
## Literatur

1. Macaluso et al. (1999) Prostate-specific antigen in vaginal fluid as a biologic marker of condom failure. *Contraception* 59: 195-201.
2. Lawson et al. (1998) Objective markers of condom failure. *Sex Transm Dis* 25: 427-423.
3. Hochmeister et al. (1999) Evaluation of Prostate-Specific Antigen (PSA) Membrane Test Assays for the Forensic Identification of Seminal Fluid: *J Forensic Sci Vol 44*: 1057-1060.
4. Hochmeister et al. (1997) Evaluation of Prostrate-Specific Antigen (PSA) Membrane Tests for the Forensic Identification of Semen. 8th International Symposium on Human Identification
5. Bolduc et al., Urinary PSA: A potential useful marker when serum PSA is between 2.5 and 10 ng/mL, *Can Urol Assoc J.* 2007 November; 1 (4): 377-381.
6. Sato et al. (2002) Use of the „SMITEST“ PSA card to identify the presence of prostate-specific antigen in semen and male urine. *Forensic Sci Int* 127: 71-74.
7. Laux et al., Forensic Detection of Semen III. Detection of PSA Using Membrane Based Tests: Sensitivity Issues with Regards to the presence of PSA in other Body Fluids.
8. Laux et al., Forensic Detection of Semen II. Comparison of the Abacus Diagnostics OneStep ABA card p30 Test and the SERATEC® PSA Semiquant Kit for the Determination of the Presence of Semen in Forensic Cases.
9. Gartside et al., Estimation of Prostate-Specific Antigen (PSA) Extraction Efficiency from Forensic Samples Using the SERATEC® PSA Semiquant Semiquantitative Membrane Test. *Forensic Science Communications* 2003 April; 5 (2).
10. SERATEC GmbH: Summary about PSA in body fluids: [http://www.seratec.com/docs/user\\_instructions/psa\\_in\\_body\\_fluids](http://www.seratec.com/docs/user_instructions/psa_in_body_fluids)
11. R. Miteva, S. Yotov, P. Georgiev, I. Fasulkov. Determination of Species Specificity of Prostate-specific Antigen in Semen. In: *Trakia Journal of Sciences* 2006, Vol. 4, No. 3, 64-68.
12. Scott McWilliams, B.S. and Bill Gartside, M.F.S. Identification of Prostate-Specific Antigen and Spermatozoa from a Mixture of Semen and Simulated Gastric Juice, *J Forensic Sci*, May 2009, Vol. 54, No. 3, 610-611.
13. Macaluso et al., Vaginal swab specimen processing methods influence performance of rapid semen detection tests: A cautionary tale *Contraception* 2010, 82(3): 291-295.

## Symbole

	Verfallsdatum
	Lagertemperatur
<b>LOT</b>	Chargennummer

Rev: 08.2013

# Instructions for use: SERATEC PSA Semiquant Cat No:PSM400F

## Intended Use

The SERATEC® PSA SEMIQUANT test is a chromatographic immunoassay for the rapid semi-quantitative detection of PSA in forensic samples. It contains two monoclonal murine anti-human-PSA antibodies as active compounds.

## Introduction

PSA (Prostate-specific antigen) is a glycoprotein produced in the prostate and secreted into the seminal fluid. PSA is one of the major proteins in seminal fluid with concentrations between 0.2 and 3.0 mg/mL. Its main function is to liquefy the seminal fluid. The high concentration of PSA in seminal fluid and its very low concentration in female vaginal fluid (0.4-0.9 ng/mL<sup>1</sup> resp. 0.0-1.25 ng/mL<sup>2</sup>, respectively) makes PSA a suitable marker in forensic casework for identifying even small amounts of seminal fluid. Particularly, these are the benefits of PSA determination in forensic biology:

- PSA can be detected in seminal fluid without spermatozoa (e.g. seminal fluid of vasectomized men).<sup>3,4</sup>
- PSA shows a high stability and could be detected in 30-year-old semen stains.<sup>3</sup>
- In vaginal smears it is detectable up to 47 hours after intercourse.<sup>1,4</sup>
- It is possible to detect PSA from vomit samples at least up to 4 h. This was shown in a study by Gartside et al. using simulated gastric juice.<sup>12</sup>
- PSA is a more specific marker than acid phosphatase.<sup>2,3</sup>
- Macaluso et al. demonstrate that the detection of PSA in vaginal swabs is more reliable than the detection of seminogelin.<sup>13</sup>

Other body fluids like blood or urine can also contain PSA, and this should be regarded when interpreting results obtained by SERATEC® PSA Semiquant. Whereas the PSA concentration in male blood serum is normally low (<4 ng/mL) and may only in the case of prostatic diseases increase up to 200 ng/mL, the PSA concentration in the urine of healthy men was reported up to 160 ng/mL with a median of 105 ng/mL.<sup>5</sup> Studies report that low amounts of PSA are already detectable in the urine of eleven years old boys.<sup>6</sup>

In case of doubt, seminal fluid and urine can be differentiated by sample dilution. PSA Semiquant will most likely react negative with urine samples diluted 1:500. Semen samples, however, generally show positive results with PSA Semiquant even diluted at 1:10<sup>6</sup>, which means a 2000 fold higher dilution. Other body fluids are containing significantly lower amounts of PSA, making an interference with PSA Semiquant results highly unlikely. A study by Laux et al. indicates that the forensic biologist can extract material from vaginal swabs and stains on clothing and be confident that a positive result is due to the presence of semen.<sup>7</sup>

**Note:** The manufacturer summarized data on the presence of PSA in body fluids along with recommendations for forensic biologists in a free available document.<sup>10</sup>

## Materials provided

- 40 individually sealed PSA Semiquant cassettes, with one plastic pipette each
- 50 mL of standard buffer solution
- user instruction leaflet

Materials required but not provided: Timer

## Storage and Stability

Both cassettes and buffer are stable up to the batch expiration dates stated on the sealed pouch and the buffer tube. Cassettes and buffer can be stored at room temperature or refrigerated (2-30 °C resp. 38-86 °F). The cassette must remain in the sealed pouch until use.

## Qualitative Characteristics

**Sensitivity:** result line is visible at as low as 1ng/mL PSA

The High dose Hook effect does not lead to false negative results: Seminal fluid is tested positive in the dilution range from 1:1 to 1:10<sup>6</sup> using the provided buffer.

## Specificity

The product shows no cross reactivity with other seminal fluid proteins. A study by Laux et al. shows that no cross reactivity was observed with neat breast milk, female urine, or female serum samples.<sup>7</sup>

No cross reactivity was observed with the seminal fluid of other mammals (dog, cat, horse, bull, pig, ram, jack)<sup>3, 11</sup> except for the seminal fluid of primates.

## Reference Preparations

The qualitative characteristics of the test are confirmed in a final QC testing using the following WHO standard: Prostate Specific Antigen NIBSC Code 96/670.

## Precautions

Forensic samples and all materials coming in contact with it should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire. The product itself and all materials coming in contact with forensic samples should be autoclaved before disposal.

- Do not use cassettes or buffer after expiration date.
- Do not use cassettes if the pouch has been damaged.
- The product consists of potentially infectious materials (e.g. antibodies). These materials do not cause any danger if the product is used according to the instructions.
- Do not open pouches until ready to perform the assay.
- Do not freeze the test cassettes.

## General Hints


- It is strongly recommended to use the provided buffer. It is adapted for obtaining optimal results. Other buffers may lead to varying line intensities.
- Do not use liquids with a pH-value below 3 or above 12 for testing. It can cause false or invalid results.
- Tissue particles do not interfere with the test result.
- Samples being subject to storage of one day or longer should be kept dry and cold (2-8 °C). Liquid samples could be frozen.
- Oral, anal and vaginal swabs shall be extracted about 2h.<sup>8</sup>
- Viscous samples should be centrifuged for about one minute at 13,000g. The PSA containing supernatant is used for testing. If the supernatant is still too viscous it should be further diluted.
- Pure semen samples and samples of very high semen concentration should be diluted at least 1:500 prior to use because of their extremely high PSA content.
- Cotton swabs and fabric or condom cuttings should be extracted in an adequate amount of the provided dilution buffer. E.g., a cotton swab or a 1 cm<sup>2</sup> cutting could be extracted in 1mL of buffer.<sup>7,9</sup>
- The addition of neat liquid samples from any source is not recommended. Liquid samples should be diluted at least 1:500 prior to testing.<sup>7</sup>

# Instructions for use: SERATEC PSA Semiquant Cat No:PSM400F

## Start of the assay

- Allow all test components to warm up to room temperature before starting the test.
- Remove the cassette out of the foil pouch and tag the cassette with a marker, if necessary.
- Add three drops (about 120 µl) in the sample well. Keep remaining sample if possible, in case it might be necessary to test additional dilutions.
- Read result after 10 minutes incubation time at room temperature. There should be no remaining fluid in the sample well at this time point.

## Interpretation of results

After **10 minutes**, PSA **negative** (below detection limit) samples will show **two lines** in the result well, whereas **PSA positive** samples will show **three lines**:


**Test result Line (T):** Reflects PSA concentration of the sample, visible in PSA positive samples only.

**Internal Standard:** Color intensity correlates with a concentration of approximately 4 ng/mL PSA.

**Control Line (C):** Control for possible procedural errors and for the integrity of test components.


**Negative result** (PSA is not detectable, either there is no PSA in the sample or the PSA concentration is below the detection limit)

Test result line (T) is not visible. Appearance of internal standard line and control line (C) confirm validity of the test.

**Positive result** (PSA detectable)


Test result line (T), internal standard line, and control line (C) appear. Any visible T-line (strong or weak colored) indicates a positive result.

## Invalid result

Internal standard line and/or control line (C) are not detectable.

In this case the assay should be repeated.

## References

1. Macaluso et al. (1999) Prostate-specific antigen in vaginal fluid as a biologic marker of condom failure. *Contraception* 59: 195-201.
2. Lawson et al. (1998) Objective markers of condom failure. *Sex Transm Dis* 25: 427-423.
3. Hochmeister et al. (1999) Evaluation of Prostate-Specific Antigen (PSA) Membrane Test Assays for the Forensic Identification of Seminal Fluid: *J Forensic Sci Vol 44*: 1057-1060.
4. Hochmeister et al. (1997) Evaluation of Prostrate-Specific Antigen (PSA) Membrane Tests for the Forensic Identification of Semen. 8th International Symposium on Human Identification at
5. Bolduc et al., Urinary PSA: A potential useful marker when serum PSA is between 2.5 and 10 ng/mL, *Can Urol Assoc J*. 2007 November; 1 (4): 377-381.
6. Sato et al. (2002) Use of the „SMITEST” PSA card to identify the presence of prostate-specific antigen in semen and male urine. *Forensic Sci Int* 127: 71-74.
7. Laux et al., Forensic Detection of Semen III. Detection of PSA Using Membrane Based Tests: Sensitivity Issues with Regards to the presence of PSA in other Body Fluids.
8. Laux et al., Forensic Detection of Semen II. Comparison of the Abacus Diagnostics OneStep ABA card p30 Test and the SERATEC® PSA Semiquant Kit for the Determination of the Presence of Semen in Forensic Cases.
9. Gartside et al., Estimation of Prostate-Specific Antigen (PSA) Extraction Efficiency from Forensic Samples Using the SERATEC® PSA Semiquant Semiquantitative Membrane Test. *Forensic Science Communications* 2003 April; 5 (2).
10. SERATEC GmbH: Summary about PSA in body fluids: [http://www.seratec.com/docs/user\\_instructions/psa\\_in\\_body\\_fluids](http://www.seratec.com/docs/user_instructions/psa_in_body_fluids)
11. R. Miteva, S. Yotov, P. Georgiev, I. Fasulkov. Determination of Species Specificity of Prostate-specific Antigen in Semen. In: *Trakia Journal of Sciences* 2006, Vol. 4, No. 3, 64-68.
12. Scott McWilliams, B.S. and Bill Gartside, M.F.S. Identification of Prostate-Specific Antigen and Spermatozoa from a Mixture of Semen and Simulated Gastric Juice, *J Forensic Sci*, May 2009, Vol. 54, No. 3, 610-611.
13. Macaluso et al., Vaginal swab specimen processing methods influence performance of rapid semen detection tests: A cautionary tale *Contraception* 2010, 82(3): 291-295.

## Symbols

	Expiration date
	Storage temperatur
<b>LOT</b>	Lot number

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1. Macaluso et al. (1999) Prostate-specific antigen in vaginal fluid as a biologic marker of condom failure. *Contraception* 59: 195-201.
2. Lawson et al. (1998) Objective markers of condom failure. *Sex Transm Dis* 25: 427-423.
3. Hochmeister et al. (1999) Evaluation of Prostate-Specific Antigen (PSA) Membrane Test Assays for the Forensic Identification of Seminal Fluid: *J Forensic Sci Vol 44*: 1057-1060.
4. Hochmeister et al. (1997) Evaluation of Prostrate-Specific Antigen (PSA) Membrane Tests for the Forensic Identification of Semen. 8th International Symposium on Human Identification
5. Bolduc et al., Urinary PSA: A potential useful marker when serum PSA is between 2.5 and 10 ng/mL, *Can Urol Assoc J*. 2007 November; 1 (4): 377-381.
6. Sato et al. (2002) Use of the „SMITEST” PSA card to identify the presence of prostate-specific antigen in semen and male urine. *Forensic Sci Int* 127: 71-74.
7. Laux et al., Forensic Detection of Semen III. Detection of PSA Using Membrane Based Tests: Sensitivity Issues with Regards to the presence of PSA in other Body Fluids.
8. Laux et al., Forensic Detection of Semen II. Comparison of the Abacus Diagnostics OneStep ABA card p30 Test and the SERATEC® PSA Semiquant Kit for the Determination of the Presence of Semen in Forensic Cases.
9. Gartside et al., Estimation of Prostate-Specific Antigen (PSA) Extraction Efficiency from Forensic Samples Using the SERATEC® PSA Semiquant Semiquantitative Membrane Test. *Forensic Science Communications* 2003 April; 5 (2).
10. SERATEC GmbH: Summary about PSA in body fluids: [http://www.seratec.com/docs/user\\_instructions/psa\\_in\\_bod\\_y\\_fluids](http://www.seratec.com/docs/user_instructions/psa_in_bod_y_fluids)
11. R. Miteva, S. Yotov, P. Georgiev, I. Fasulkov. Determination of Species Specificity of Prostate-specific Antigen in Semen. In: *Trakia Journal of Sciences* 2006, Vol. 4, No. 3, 64-68.
12. Scott McWilliams, B.S. and Bill Gartside, M.F.S. Identification of Prostate-Specific Antigen and Spermatozoa from a Mixture of Semen and Simulated Gastric Juice, *J Forensic Sci*, May 2009, Vol. 54, No. 3, 610-611.
13. Macaluso et al., Vaginal swab specimen processing methods influence performance of rapid semen detection tests: A cautionary tale *Contraception* 2010, 82(3): 291-295.

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1. Macaluso et al. (1999) Prostate-specific antigen in vaginal fluid as a biologic marker of condom failure. *Contraception* 59: 195-201.
2. Lawson et al. (1998) Objective markers of condom failure. *Sex Transm Dis* 25: 427-423.
3. Hochmeister et al. (1999) Evaluation of Prostate-Specific Antigen (PSA) Membrane Test Assays for the Forensic Identification of Seminal Fluid: *J Forensic Sci Vol 44: 1057-1060.*
4. Hochmeister et al. (1997) Evaluation of Prostrate-Specific Antigen (PSA) Membrane Tests for the Forensic Identification of Semen. 8th International Symposium on Human Identification at
5. Bolduc et al., Urinary PSA: A potential useful marker when serum PSA is between 2.5 and 10 ng/mL, *Can Urol Assoc J.* 2007 November; 1 (4): 377-381.
6. Sato et al. (2002) Use of the „SMITEST” PSA card to identify the presence of prostate-specific antigen in semen and male urine. *Forensic Sci Int* 127: 71-74.
7. Laux et al., Forensic Detection of Semen III. Detection of PSA Using Membrane Based Tests: Sensitivity Issues with Regards to the presence of PSA in other Body Fluids.
8. Laux et al., Forensic Detection of Semen II. Comparison of the Abacus Diagnostics OneStep ABA card p30 Test and the SERATEC® PSA Semiquant Kit for the Determination of the Presence of Semen in Forensic Cases.
9. Gartside et al., Estimation of Prostate-Specific Antigen (PSA) Extraction Efficiency from Forensic Samples Using the SERATEC® PSA Semiquant Semiquantitative Membrane Test. *Forensic Science Communications* 2003 April; 5 (2).
10. SERATEC GmbH: Summary about PSA in body fluids: [http://www.seratec.com/docs/user\\_instructions/psa\\_in\\_body\\_fluids](http://www.seratec.com/docs/user_instructions/psa_in_body_fluids)
11. R. Miteva, S. Yotov, P. Georgiev, I. Fasulkov. Determination of Species Specificity of Prostate-specific Antigen in Semen. In: *Trakia Journal of Sciences* 2006, Vol. 4, No. 3, 64-68.
12. Scott McWilliams, B.S. and Bill Gartside, M.F.S. Identification of Prostate-Specific Antigen and Spermatozoa from a Mixture of Semen and Simulated Gastric Juice, *J Forensic Sci*, May 2009, Vol. 54, No. 3, 610-611.
13. Macaluso et al., Vaginal swab specimen processing methods influence performance of rapid semen detection tests: A cautionary tale *Contraception* 2010, 82(3): 291-295.

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1. Macaluso et al. (1999) Prostate-specific antigen in vaginal fluid as a biologic marker of condom failure. *Contraception* 59: 195-201.
2. Lawson et al. (1998) Objective markers of condom failure. *Sex Transm Dis* 25: 427-423.
3. Hochmeister et al. (1999) Evaluation of Prostate-Specific Antigen (PSA) Membrane Test Assays for the Forensic Identification of Seminal Fluid: *J Forensic Sci Vol* 44: 1057-1060.
4. Hochmeister et al. (1997) Evaluation of Prostrate-Specific Antigen (PSA) Membrane Tests for the Forensic Identification of Semen. 8th International Symposium on Human Identification
5. Bolduc et al., Urinary PSA: A potential useful marker when serum PSA is between 2.5 and 10 ng/mL, *Can Urol Assoc J.* 2007 November; 1 (4): 377-381.
6. Sato et al. (2002) Use of the „SMITEST“ PSA card to identify the presence of prostate-specific antigen in semen and male urine. *Forensic Sci Int* 127: 71-74.
7. Laux et al., Forensic Detection of Semen III. Detection of PSA Using Membrane Based Tests: Sensitivity Issues with Regards to the presence of PSA in other Body Fluids.
8. Laux et al., Forensic Detection of Semen II. Comparison of the Abacus Diagnostics OneStep ABA card p30 Test and the SERATEC® PSA Semiquant Kit for the Determination of the Presence of Semen in Forensic Cases.
9. Gartside et al., Estimation of Prostate-Specific Antigen (PSA) Extraction Efficiency from Forensic Samples Using the SERATEC® PSA Semiquant Semiquantitative Membrane Test. *Forensic Science Communications* 2003 April; 5 (2).
10. SERATEC GmbH: Summary about PSA in body fluids: [http://www.seratec.com/docs/user\\_instructions/psa\\_in\\_body\\_fluid\\_s](http://www.seratec.com/docs/user_instructions/psa_in_body_fluid_s)
11. R. Miteva, S. Yotov, P. Georgiev, I. Fasulkov. Determination of Species Specificity of Prostate-specific Antigen in Semen. In: *Trakia Journal of Sciences* 2006, Vol. 4, No. 3, 64-68.
12. Scott McWilliams, B.S. and Bill Gartside, M.F.S. Identification of Prostate-Specific Antigen and Spermatozoa from a Mixture of Semen and Simulated Gastric Juice, *J Forensic Sci*, May 2009, Vol. 54, No. 3, 610-611.
13. Macaluso et al., Vaginal swab specimen processing methods influence performance of rapid semen detection tests: A cautionary tale *Contraception* 2010, 82(3): 291-295.

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