

SERATEC Product List 2010

Product Line	Product Name	Cat.No.	Specimen	Format	Width	Cut-Off	Standard Kit Size*	Application	CE
Pregnancy & Fertility	hCG Cassette Urine	HCG-C1	Urine	Cassette	4,00 mm	25 mIU/mL	Box, 40 tests	professional	Yes
	hCG Cassette Combo	HCG-C2	Urine / Serum	Cassette	4,00 mm	25 mIU/mL	Box, 40 tests	professional	Yes
	hCG Strip, single sealed	HCG-S1	Urine	Strip	5,00 mm	25 mIU/mL	Box, 25 tests	professional	Yes
	hCG Strip in tube	HCG-S2	Urine	Strip	5,00 mm	25 mIU/mL	Tube, 20 tests	professional	Yes
	hCG Strip small, in tube	HCG-S3	Urine	Strip	2,50 mm	25 mIU/mL	Tube, 50 tests	professional	Yes
	hCG Direct	HCG-D1	Urine	Direct	4,00 mm	25 mIU/mL	Packet, 1 test	home	Yes
	hCG Midstream	HCG-M1	Urine	Midstream	8,00 mm	25 mIU/mL	Packet, 1 test	home	Yes
	hCG Midstream, early hCG test	HCG-M2	Urine	Midstream	8,00 mm	15 mIU/mL	Packet, 1 test	home	Yes
	hCG Strip, sealed in envelope	HCG-S4	Urine	Strip	8,00 mm	25 mIU/mL	Envelope, 1 test	home	Yes
	hCG Strip, early hCG test	HCG-S5	Urine	Strip	5,00 mm	15 mIU/mL	Envelope, 1 test	home	Yes
	hCG Strip, early & confirmative	HCG-S45	Urine	Strip	5,00 mm	15 & 25 mIU/mL	Envelope, 2 tests	home	Yes
	LH Max Ovulation Cassette	LH-C1	Urine	Cassette	8,00 mm	35 ng/mL	Box, 6 tests	home	Yes
	LH Max Ovulation Strip	LH-S1	Urine	Strip	5,00 mm	35 ng/mL	bulk	home	Yes
LH Max Ovulation Midstream	LH-M1	Urine	Strip	5,00 mm	35 ng/mL	bulk	home	Yes	
Drugs of Abuse	Drug Screen Amphetamine	DSC-AMP	Urine	Cassette	4,00 mm	500 / 1000 ng/mL	Box, 40 tests	professional	Yes
	DSS-AMP	Urine	Strip	4,00 mm	1000 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Barbiturates	DSC-BAR	Urine	Cassette	4,00 mm	300 ng/mL	Box, 40 tests	professional	Yes
	DSS-BAR	Urine	Strip	4,00 mm	300 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Benzodiazepine	DSC-BZO	Urine	Cassette	4,00 mm	300 ng/mL	Box, 40 tests	professional	Yes
	DSS-BZO	Urine	Strip	4,00 mm	300 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Buprenorphine	DSC-BUP	Urine	Cassette	4,00 mm	10 ng/mL	Box, 40 tests	professional	Yes
	DSS-BUP	Urine	Strip	4,00 mm	10 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Cocaine	DSC-COC	Urine	Cassette	4,00 mm	300 ng/mL	Box, 40 tests	professional	Yes
	DSS-COC	Urine	Strip	4,00 mm	300 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen MDMA	DSC-MDMA	Urine	Cassette	4,00 mm	500 ng/mL	Box, 40 tests	professional	Yes
	DSS-MDMA	Urine	Strip	4,00 mm	500 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Methadone	DSC-MTD	Urine	Cassette	4,00 mm	300,00 mm	Box, 40 tests	professional	Yes
	DSS-MTD	Urine	Strip	4,00 mm	300,00 mm	Box, 25 tests	professional	Yes	
	Drug Screen Morphine	DSC-MOR	Urine	Cassette	4,00 mm	300 ng/mL	Box, 40 tests	professional	Yes
	DSS-MOR	Urine	Strip	4,00 mm	300 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Metamphetamine	DSC-MET	Urine	Cassette	4,00 mm	1000 ng/mL	Box, 40 tests	professional	Yes
	DSS-MET	Urine	Strip	4,00 mm	1000 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Phencyclidine	DSC-PCP	Urine	Cassette	4,00 mm	25 ng/mL	Box, 40 tests	professional	Yes
	DSS-PCP	Urine	Strip	4,00 mm	25 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Tricyclic Antidepressants	DSC-TCA	Urine	Cassette	4,00 mm	1000 ng/mL	Box, 40 tests	professional	Yes
	DSS-TCA	Urine	Strip	4,00 mm	1000 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Marijuana (THC)	DSC-THC	Urine	Cassette	4,00 mm	50 ng/mL	Box, 40 tests	professional	Yes
	DSS-THC	Urine	Strip	4,00 mm	50 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Cotinine	DSC-COT	Urine	Cassette	4,00 mm	200 ng/mL	Box, 40 tests	professional	Yes
	DSS-COT	Urine	Strip	4,00 mm	200 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen EDDP	DSC-EDDP	Urine	Cassette	4,00 mm	100 ng/mL	Box, 40 tests	professional	Yes
	DSS-EDDP	Urine	Strip	4,00 mm	100 ng/mL	Box, 25 tests	professional	Yes	
	Drug Screen Fentanyl	DSC-FEN	Urine	Cassette	4,00 mm	10 ng/mL	Box, 40 tests	professional	Yes
	DSS-FEN	Urine	Strip	4,00 mm	10 ng/mL	Box, 25 tests	professional	Yes	
Drug Screen Ketamine	DSC-KET	Urine	Cassette	4,00 mm	1000 ng/mL	Box, 40 tests	professional	Yes	
DSS-KET	Urine	Strip	4,00 mm	1000 ng/mL	Box, 25 tests	professional	Yes		
Drug Screen Tramadol	DSC-TRA	Urine	Cassette	4,00 mm	100 ng/mL	Box, 40 tests	professional	Yes	
DSS-TRA	Urine	Strip	4,00 mm	100 ng/mL	Box, 25 tests	professional	Yes		
Drug Screen Propoxyphen	DSC-PRO	Urine	Cassette	4,00 mm	300 ng/mL	Box, 40 tests	professional	Yes	
DSS-PRO	Urine	Strip	4,00 mm	300 ng/mL	Box, 25 tests	professional	Yes		
Drug Screen Oxycodone	DSC-OXY	Urine	Cassette	4,00 mm	100 ng/mL	Box, 40 tests	professional	Yes	
DSS-OXY	Urine	Strip	4,00 mm	100 ng/mL	Box, 25 tests	professional	Yes		
Drug Screen Methaqualone	DSC-MAQ	Urine	Cassette	4,00 mm	300 ng/mL	Box, 40 tests	professional	Yes	
DSS-MAQ	Urine	Strip	4,00 mm	300 ng/mL	Box, 25 tests	professional	Yes		
Drug Screen Multitest	DSMC-...	Urine	Multi-Window	-	-	Box, 40 tests	professional	Yes	
(various parameter combinations)	DSMD-...	Urine	Multi-Diptest	-	-	Box, 40 tests	professional	Yes	
Forensic Biology	PSA Semiquant	PSM400F	Traces	Cassette	5,00 mm	2 ng/mL	Box, 40 tests	professional	N/A**
	PSA Semiquant Crime Scene Kit	PSM400FE	Traces	Cassette	5,00 mm	2 ng/mL	Box, 20 tests	professional	N/A**
	HemDirect Hemoglobin Test	HBFO7	Traces	Cassette	5,00 mm	15 ng/mL	Box, 40 tests	professional	N/A**
Other	PSA Semiquant Tumor Marker	PSM400T	Serum	Cassette	5,00 mm	2 ng/mL	Box, 40 tests	professional	N/A
	PSA Antibody 28/A4	AB28-A4	-	-	-	-	-	professional	N/A
	PSA Antibody 43/A2	AB43-A2	-	-	-	-	-	professional	N/A
	SERAQUANT	SQNT	-	-	-	-	-	professional	N/A**

*The kit size, design etc. can be subject to change due to customer's requirements

** For forensic use no CE registration necessary