


Verfügbare Parameter	Nachweisgrenzen (ng/ml)		
	1000	500	300
Amphetamin	1000	500	300
Barbiturate	300	200	
Benzodiazepine	300	200	100
Buprenorphin	10		
Cotinin	600	300	200
EDDP	100		
Fentanyl	200		
Ketamin	1000		
Kokain	300	200	100
MDMA	1000	500	
Metamphetamin	1000	300	
Methadon	300		
Morphin	300	200	100
Opiate	2000		
Oxycodon	300	100	
Phencyclidin	25		
Propoxyphen	300		
THC	200	150	50
Tramadol	300	100	
Trizyklische Antidepressiva	1000		



SERATEC® Gesellschaft für Biotechnologie mbH
 Postfach 3706, D-37027 Göttingen

Tel.: +49 (0)551 504800 E-Mail: contact@seratec.com
 Fax: +49 (0)551 5048080 Web: www.seratec.com