

- No further hardware needed
- Immediate result, for both single and multi-tests up to 7 parameters
- Objectivity and result accuracy
- Intuitive and easy use
- Geotagging
- Automated report generation
- PDF export
- XML export enables integration with lab management software


»Intuitive & accurate scanning
and reporting solution«



Serascope

Parameters Available Standard Cut-Off (ng/mL)

Amphetamine	1000	500	300
Barbiturates	300	200	
Benzodiazepines	300	200	100
Buprenorphine	10		
Cocaine	300	200	100
Cotinine	600	300	200
EDDP	100		
Fentanyl	200		
Ketamine	1000		
MDMA	1000	500	
Metamphetamine	1000	300	
Methadone	300		
Morphine	300	200	100
Opiates	2000		
Oxycodone	300	100	
Phencyclidine	25		
Propoxyphene	300		
THC	200	150	50
Tramadol	300	100	
Tricyclic Antidepressants	1000		



SERATEC® Gesellschaft für Biotechnologie mbH
P.O. Box 3706, D-37027 Goettingen, Germany

Phone: +49 (0)551 504800 E-Mail: contact@seratec.com
Fax: +49 (0)551 5048080 Web: www.seratec.com